



A division of Gene Therapy Systems, Inc.

siGuard™ RNase Inhibitor

PRODUCT SUMMARY

Cat. No: T520100

Description: siGuard RNase Inhibitor provides optimal protection of siRNA from degradation due to contamination by ribonucleases.

FOR RESEARCH USE ONLY

Components:

siGuard RNase Inhibitor (1 vial x 100 units in 100 µl)

Shipping and Storage: The siGuard RNase Inhibitor is shipped frozen. For maximum stability and long-term use, store at -20°C upon receipt. The siGuard RNase Inhibitor is stable for six months when stored properly.

Unit Definition: one unit is the amount of RNase inhibitor required to inhibit by 50% the activity of 5 ng of ribonuclease A at 25°C. Assay Conditions: 0.1M Tris-acetate, pH 6.5, 1mM EDTA, 1.0 mM cyclic 2',3'-cytidine monophosphate. Reaction volume 1.0 ml. Inhibitor activity is determined by the inhibition of RNase A hydrolysis of cyclic 2',3'-cytidine monophosphate in a spectrophotometric recording assay.

Storage Buffer: 20 mM HEPES-KOH, pH 7.6, 50 mM KCl, 5mM dithiothreitol and 50% glycerol.

Molecular Weight: ~50,000 by SDS-PAGE.

Optimal pH: 5-9 (highest activity at pH 7-8).

METHODS AND PROCEDURES

1. Use 1 unit for 250 ul of transfection volume. Please refer to Table 1 for example transfection reactions.
2. Denaturing conditions (*i.e.* urea, temperatures > 50°C) should be avoided as they may cause release of active ribonuclease from the complex.

3. Ribonuclease inhibitor can be removed by phenol extraction or inactivated by heating at 65°C for 10 minutes.

Table 1: Suggested siGuard Volumes for siRNA Transfection

Tissue culture plate size	Total transfection volume / well	Amount of siGuard / well
96 wells	100 µl	0.5 µl
48 wells	200 µl	1 µl
24 wells	500 µl	2 µl
6 wells	1000 µl	4 µl

References:

Blackburn, P. (1979) Ribonuclease inhibitor from human placenta: rapid purification and assay *J.Biol.Chem.* **254** (24): 12484-12487.

RELATED PRODUCTS

Product	Cat. No.
<i>For efficient generation of siRNA</i>	
Dicer siRNA Construction Kit 5 genes / 50 transfections	T510001
<i>For efficient and functional siRNA transfection</i>	
GeneSilencer® siRNA Transfection Reagent 200 reactions (0.75 ml)	T500750
GeneSilencer® siRNA Transfection Reagent 5 x 200 reactions (5 x 0.75 ml)	T505750

Please contact us for our latest catalogue of life science research products.

Genlantis

A division of Gene Therapy Systems, Inc.
10190 Telesis Court
San Diego, CA 92121
Tel: 858-457-1919 • 888-428-0558
Fax: 858-623-9494

Website: <http://www.genlantis.com>