
20000004261

S 65 KG_032015

IKA®

Dispergierwerkzeug S 65 KG-HH-G 65 G/M/F
Dispersion tool S 65 KG-HH-G 65 G/M/F
Outil dispersant S 65 KG-HH-G 65 G/M/F
Диспергирующий инструмент S 65 KG-HH-G 65 G/M/F

Betriebsanleitung	DE	3
Operating instructions	EN	7
Mode d'emploi	FR	11
Инструкция по эксплуатации	RU	15


Fig. 1

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Content

	Page
Explication of warning symbols	7
Safety instructions	7
Product information	7
Application instruction	8
Assembly dispersing tools	8
Cleaning, disinfection and sterilization	8
Maintenance	9
Spare parts	10
Technical Data	10
Warranty	10

Explication of warning symbols

General hazard

This symbol identifies information **that is of vital importance for protecting your health and safety**. Disregarding this information may lead to health impairment and injuries.

This symbol identifies information **that is of importance for the technically correct functioning of the system**. Disregarding this information may result in damage to the device or to system components.

Safety instructions

• **Read the operating instructions in full before starting up and follow the safety instructions.**

- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. Otherwise there is a risk of splashing and evaporation of liquids; body parts, hair, clothing and jewellery getting caught.
- Use the dispersing tool always inside the stirring vessel.
- Never run dispersion tools dry, as the gasket and bearings will be destroyed if the tools are not cooled by the medium.
- Rapid temperature change (thermal shock!) can destroy the hard metal sealing surfaces of the dispersing tool **S 65-KG HH**.
- No liquid may get into the drive unit of the disperser.
- The work has to be interrupted immediately if you notice unusual noise. The bearings and the seals must be controlled, and replaced if necessary.
- Note the operating instructions of the disperser unit.
- Abrasion of the dispersion equipment or the rotating accessories can get into the medium you are working on.
- If a dispersing element is damaged, e.g. with hair-line cracks, chips or bent rotor/stator teeth, it must not be used.

Product information

- The dispersing tool consists of a shaft and a generator (stator/rotor). The shaft carries the marking **S 65 KG-HH**, where **KG** stands for the ball bearing and **HH** for a slide ring seal made of hard metal. The seal prevents liquid rising up the shaft and destroying the ball bearing. The top end of the shaft is in the form of a plug-in coupling. The bottom end of the shaft carries a thread for mounting the stator. The shaft – running in ball bearings- carries the rotor.
- **Content of package:** Shaft, dispersion head and operating instructions.
- For the **S 65** range, there are three types of generator, designated **G 65 G**, **G 65 M** and **G 65 F**, where the prefix **G** stands for generator and the suffices **G** for coarse, **M** for medium and **F** for fine. The appropriate generator must be used according to the medium. If the feed material is too coarse the generator can become clogged. The sealing faces then run dry and can be damaged. The generators are suitable for solids components with the following grain sizes:

Generator type	Grain size
G 65 G	40 mm
G 65 M	8 mm
G 65 F	1 mm

It can often be helpful to process a medium with all three generators in sequence, starting with generator **G 65 G**. In such circumstances, to avoid having to replace the entire generator it is worthwhile procuring three shafts, which can be replaced using just a few action steps.

Application instruction

- Remove the protection cap (1) before use. Use this cap to protect shaft tubes from damage when not in service..
- **Attention!** Never run the dispersion tool dry. The sealing surfaces will be destroyed without cooling by the medium.
- The dispersion tools are not suitable for continuous operation.
- Dispensing causes the medium to heat up.
- Likewise the optimal dispersion duration and rotating frequency must be determined by attempts.
- The function of the dispersion tools depends on the condition of the sharp edges of rotor and stator. Abrasive media can round off these edges fast, whereby the disperse effect diminishes.
- **Application:** Wet crushing, dispersions, emulsifying (Batch-operation).

Assembly dispersing tools

Pre-assembly:

Before the dispersion tools are connected to the dispersion unit, the shaft and generator are pre-assembled. The following tools (set of tools included in the **T 65 basic/digital ULTRA-TURRAX®** delivery package) are required to change the generator:

Rotor wrench

Hook wrench

Shaft wrench


Right-angle screwdriver DIN 911 AF 4


The slide ring seal must be correctly fitted before pre-assembly of the shaft and generator!


Assemble the generator on the shaft as follows:

Screw the stator (30) on to the stator tube (15) counterclockwise by hand. Note the left-hand thread.

Gently tighten the stator (30) with the hook wrench. The stator is turned in the same direction later, by the flow pressure of the rotor (40), and hence cannot come loose.

Screw the rotor (40) clockwise on to the shaft (16) (right-hand thread). To do this, hold the shaft (16) with the shaft wrench at the coupling section and tighten the rotor (40) with the rotor wrench. With the **G 65 G** generator, the rotor wrench is inserted between the teeth of the stator and turned.

Mounting the dispersing tool on the unit:

The dispersion tool must only be mounted on the dispersion unit when the mains plug is pulled out and the drive shaft is stationary.

Insert the dispersing tool in the dispersion tool receptacle and secure it with the handwheel bolt. Make sure that the shaft is inserted right up to the stop in the receptacle. It can be helpful to turn the shaft lightly as it is inserted. The tool is locked in position by tightening the handwheel bolt. Make sure that the handwheel bolt is securely tightened, and also check the tightness frequently when the unit is in service.

Cleaning, disinfection and sterilization


Cleaning the shaft and generator:

The shaft must be cleaned immediately after the unit has been used. This prevents any substance residues adhering to the rotor and stator threads, or bacteria cultures forming in undesired places. Run the dispersion tool in a solvent which will dissolve substance residues but not harm the seals. Because the flow speeds are so high, in such circumstances the rotor and stator are quickly cleaned.

Do not put the dispersion tool in a vessel containing solvent to dissolve substance residues. The seals are mainly dynamic-acting, and quickly start to leak when stationary.

Clean the dispersion tools and shaft as follows:

1. With the drive switched off, dismount the dispersion tool from the drive by slackening the handwheel bolt.
2. Unscrew the generator from the shaft. The disassembly procedure is the reverse of the assembly procedure (see “**Pre-assembly of dispersion tool**”).

The stator can be very tightly seated on the shaft tube, particularly when working with viscous substances, or when cleaning has been neglected, or when medium has penetrated into the thread between shaft tube and stator. In this case, the shaft tube can be held in a vice with soft jaws to crack the thread. Note the left-hand thread!

3. Rinse all threads and clean them with a brush if necessary.

The friction surfaces of the seal must be clean.

The shaft and dispersion tool can now be preassembled again, ready for use. Reassemble the dispersion tool in the reverse sequence. Attention! Do not tighten the rotor too strongly, because otherwise the teeth at the rotor become bent and the dispersion tool will be destroyed.

Sterilization and disinfection:

Wet chemical methods are permissible for sterilization and disinfection of the dispersion tools. Many instances of disinfection can be solved by means of bactericidal solutions. It is important that the remains of the disinfectant are subsequently removed with sterilized water.

No moisture must be allowed to enter the ball bearing, so, basically, the upper end of the shaft tube must be sealed (e.g. with silicone plugs).

Other sterilisation processes (such as steam sterilisation or disinfection by hot air sterilisation) are not permitted.

Maintenance

The shaft seals must be regularly checked by dismantling them and examining the contact faces to ensure they are not cracked or pitted. If there is a leak at the seal the suction effect can cause the medium to rise up in the tube around the rotating shaft and may penetrate the drive unit. If the contact faces are cracked or pitted, the rotary seal must be replaced.

Removal and fitting of the rotary seal (see Fig. 1):

To prevent damage to the unit, the slide ring seal must only be changed by a skilled mechanical person.

Change the slide ring seal as follows:

1. Using a 4 AF hexagon L-wrench to DIN 911, undo the three socket-head threaded pins (3). The lower part of the rotary seal (18b) can then easily be withdrawn.
2. To remove the drive shaft, the upper circlip (8) must be removed using a pair of circlip pliers. Now the drive shaft can be pushed upwards and out the shaft tube. After this, the upper part of the rotary seal (18a) can be extracted from the shaft tube.
3. When refitting, adjust the bearing pressure of the two sealing surfaces. The correct contact pressure is obtained when the inner part of the lower half of the rotary seal (18c) is pushed in sufficiently far before tightening that it is flush with the outer part of the lower half of the rotary seal (18b).

Spare parts

Item	Designation
3	Notched knurled pin
4	Grooved ball bearing
5	Grooved ball bearing
6	Circlip
8	Circlip
15	Shaft tube
16	Shaft
18	Rotary seal
2001	Rotary seal, compl.
19	O-ring
20	O-ring
21	O-ring
30	Stator
40	Rotor

Spare parts diagram see front page:

For orders of spare parts please specify the fabrication number on the shaft tube, the type of dispersion tool and the designation of the spare part.

Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard.

For this use the “**Decontamination Certificate**” included or download the form from the **IKA®** website **www.ika.com**.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Technical Data

Designation	S65KG-HH-G 65G	S65KG-HH-G65M	S65KG-HH-G65F
Stator/Rotor (Ø)	65/58 mm		
Max. permissible speed	9500 rpm		
Working range pressure	6 bar		
Working range vacuum	1mbar		
Peripheral speed at 7200 rpm	21.9 m/s		
Min./max. immersion depth	90 - 450 mm	80 - 450 mm	80 - 450 mm
Material in contact with medium	FFPM / SIC, AISI 316L		
Suitable for solvents	yes		
Temperature range upto	180 °C		
Ultimate fineness dispersions	25 - 75 µm	20 - 50 µm	5 - 20 µm
Ultimate fineness emulsions	5 - 25 µm	5 - 15 µm	1 - 10 µm

Warranty

In accordance with **IKA®** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

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IKA®-Werke GmbH & Co.KG