General instruction for STCON7 electrode

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Specifications:

Connection:	Min-Din
Cable length:	1.0m
Shaft length:	95mm
Shaft diameter:	12mm
Shaft material:	ABS
Cell material:	316L steel
Temperature range:	0-60°C
Measurement range:	0.02~200 μ S/cm
Accuracy:	0.02 μ S/cm

This STCON7 Conductivity Measuring Cell utilizes the 2-pole potentiometric method for measuring conductivity. This electrode is especially designed for low conductivity (pure water, distilled water etc.) measurement in samples.

Furthermore, the combined temperature sensor makes STCON7 automatic temperature compensating.

It is a low conductivity probe, that is ideal for lab use.

Operating Instructions

For optimum performance, use the following operating procedures:

1. Remove all packaging material from the probe before use. Retain the quality certificate and the operating instructions for reference.
2. To prevent carry-over from solution to solution, rinse the probe with distilled water between measurements.
3. When measuring, make sure the solution is above the cell chamber.
4. Ensure the cell chamber is free of bubbles when measuring. To reduce air bubbles, immerse the probe at an angle and then rise to vertical position. Remaining small bubbles can be removed by gently tapping against the probe shaft.
5. Allow sufficient time for the probe to stabilize when measuring samples at different temperatures.
6. Clean the probe with distilled water after use and keep it dry when not in use.
7. If there is a build-up of solids inside the cell chamber, carefully remove it with a cotton bud soaked in detergent solution, and then rinse the probe with distilled water.

8. Ensure the cable connectors are free of moisture and dirt.
9. For calibration and measurement procedures, refer to the Instruction manual supplied with conductivity meter.
10. In general, measurements using stirred samples are preferred.

Storage

Rinse or clean the electrode after each use. Then keep it in a dry clean environment.

Precautions and Limitations

1. Do not expose the shaft to organic solvents, when cleaning or when taking measurements.
2. Do not use the probe outside the recommended temperature range.
3. Please calibrate the electrode with 84 $\mu\text{S}/\text{cm}$ standard solution for an accurate measurement according to the meter instruction manual.
4. STCON7 can work with Starter 300C and Starter 3100C meters.

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Ohaus Corporation, all rights reserved/ todos