

PARI SINUS™ Configurations


Configuration	Part #	NDC #
PARI SINUS™ Pulsating Aerosol System	#028F1005	#83490-0280-08
PARI LC® Sprint Sinus Nebulizer Kit	#028F2861	#83490-0280-11

Features


- Targeted delivery of medication directly to the sinuses
- Deposits significantly more medication into the sinus cavities vs. nasal sprays<sup>1</sup>
- Less medication required = lower side effects

The Difference is in the Delivery

Aerosol Delivery to the Sinuses



Non-Pulsating Aerosol



Pulsating Aerosol



Performance with LC® Sprint Sinus<sup>2</sup>

Total Output Rate = **220** mg/min  
MMD (mass median diameter) = **3.2** µm  
Mass Percentage Below 5µm = **71**%  
Maximum Fill Volume = **8** mL

Accessories:

Filter Replacement 3 Pack = **028F0012P3**  
Soft Nose Plug 3 Pack = **028F0576**

PARI SINUS

Pulsating Aerosol System



Good Vibrations.

The PARI SINUS™ Pulsating Aerosol System combines efficient nebulization with a vibrating pulse, delivering aerosol directly to the sinuses.

*From the leader in clinically proven aerosol delivery*

<sup>1</sup>Schuschnig et al. (2006) "Comparison of delivery efficiency in a nasal cast model of fluticasone (Flutide®) nasal spray versus a novel solution aerosolized via the PARI SINUS™". OTO Conference Poster 2006  
<sup>2</sup>Measured with the Malvern Spraytec at 23°C, 50% relative humidity: 0.9% NaCl (5ml), Jet Flow 4,6 l/min).