

Veterinary Pulse Oximeter

Oxy9Wave Vet



High Performance

- * 3.2" Color TFT LCD for clear view
- * Low perfusion SpO₂ measurement
- * Nellcor compatibility
- * Li-ion rechargeable battery allows cordless use for up to 20 hours

User Convenience

- * Landscape & portrait views
- * Intuitive operation
- * SpO₂, BPM, pulse parameter display



Portrait display



Landscape display

SPECIFICATIONS

■ Performance

SpO₂ (Blood oxygen saturation)	Range : 40 to 100% Accuracy : 70~100%, ±2% 0~69% Unspecified
BPM (Heartbeat)	Range : 30 to 250 bpm Accuracy : ±2bpm
Trend	Memory : 30 days at 10 second intervals

■ Display

Type	3.2" TFT LCD (Resolution : 240 × 400 pixel)
Display parameter	SpO ₂ , BPM, SpO ₂ limit, BPM limit, Level Bar, Heartbeat wave, Time and Date, Alarm volume, Battery status and brightness
LED display	Sensor separation, alarm silence, charging status and power connection

■ Power Supply

Power consumption	24VA (max)
Adaptor	Input : 100-240V~, 50-60Hz, 0.4A Output : +5V, 2.0A
Battery	Type : Lithium-ion battery (4400mAh, 3.7V) Operating time : Up to 20 hours Charging time : 6.5 hours

■ General

Operating temperature	5°C to 40°C
Non-operating temperature	-20°C to 70°C
Humidity	5% to 95%
Operating altitude	-1,000 ft to 12,000 ft
Weight	300g
Size	45(H) × 210(W) × 75(D)mm

■ Accessory

Standard	Reusable multisite SpO ₂ probe AC adaptor Pole mount bracket	1ea 1ea 1ea
Option	Transflectance SpO ₂ probe Extension cable	



South Korea: Bionet Co., Ltd. www.ebionet.com
Tel: +82-70-7585-6418 (Int'l)/+82-2-6300-6483 (Domestic)
E-mail: sales@ebionet.com (Int'l)/yunms@ebionet.com (Domestic)

North America (U.S.A): Bionet America, Inc.
Tel: +1-714-734-1760
E-mail: sales@BionetUS.com www.BionetUS.com

China : Bionet China, Inc. (中国佰奈特)
Tel: +86-10-8476-3601
E-mail: ayw1969@ebionet.com www.ebionet.com.cn

Latin America: Bionet Latin America, Inc.
Tel: +1-949-800-9585
E-mail: jayshim@ebionet.com www.ebionet.com

Europe: MGB Endoskopische Geräte GmbH Berlin (Germany)
Tel: +49-30-6392-7000
E-mail: sales@mgb-berlin.de www.mgb-berlin.de